



DEVICE CONNECTIVITY FOR  
POINT-OF-CARE TESTING

# ID NOW™ CONNECTIVITY

The ID NOW™ system is a rapid, instrument-based, isothermal system for the qualitative detection of infectious diseases. When interfaced to a RALS™ system, it provides users with numerous benefits.

## CONNECTIVITY TO RALS ENABLES USERS TO:


- Centralize all POC test results into the RALS system
- Transfer results to an EMR/EHR/LIS
- Email results to operators and managers
- Access and manage mis-identified patient and out of range QC results
- Monitor and manage operator certifications



VISIT [RALS.COM](http://RALS.COM)

# RALS™ | ID NOW™ Connectivity

## HOSPITAL / FACILITY / INSTITUTION NETWORK



### DATA FLOW

The ID NOW transfers data to the RALS System, where it is evaluated and either sent directly to an EMR/EHR/LIS or held for review.

### SUPPORTED FEATURES

- Results Management
- Operator Management
- Certification Management
- Device Management
- Patient, QC, Proficiency Management
- Evaluation Criteria Settings

### HOW THE DEVICE CONNECTS

- Direct network connection

### BI-DIRECTIONAL COMMUNICATION

Data sent from ID NOW to RALS

- Patient results
- QC results
- Test materials
- Patient ID
- Operator ID
- Date/time (performed and downloaded)
- Location

Data sent from RALS to ID NOW

- Operator List (up to 1,000 operators)
- Date and Time
- Patient Query Support\*:
  - Name
  - Date of Birth
  - Query Status (SUCCESS, NOT FOUND, MULTIPLE)
- Cartridge Material List Upload\*
- Device Configuration Support\*
  - Meter Settings
    - Positive Patient ID Enabled/Disabled
  - QC Schedule Settings

\* Available in RALS Version 6.0, SP5 and ID NOW instruments running firmware version 6.1.22

### AVAILABLE REPORTS

- Anomalous Results by Device
- Anomalous Results by Operator
- Operator Statistics
- Operator Pending Recertification Range
- Operator Pending Recertification Date
- Operator with No Certifications
- QC Test Results
- Patient Results by Operator
- Patient Results by Patient
- LIS Statistics by Location
- LIS Statistics by Operator
- Institution Statistics
- Location Statistics
- RALS System Setup
- RALS Device Log
- Workload Summary
- Test Materials

TO LEARN MORE, VISIT [RALS.COM](https://www.rals.com)

© 2020 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners. Any photos displayed are for illustrative purposes only. Any person depicted in such photos is a model. AI-COLL-0032 Rev. 6 07/20

