



DEVICE CONNECTIVITY FOR  
POINT-OF-CARE TESTING

# AFINION™ CONNECTIVITY

Afinion™ Analyzers are designed to enable quick and easy on-the-spot testing of HbA1c and Albumin/Creatinine Ratio (ACR) urine tests. When connected to a RALS™ system, users in hospitals and physician office labs can access numerous benefits.

*Connectivity by the RALS System enables users to:*


- Centralize all test results into the RALS system
- Transfer results to an EMR/EHR/LIS
- Email results to operators and managers
- Access and manage mis-identified patient and out of range QC results
- Monitor and manage operator certifications



VISIT [RALS.COM](http://RALS.COM)

# RALS™ | Afinion™ Connectivity

## HOSPITAL / FACILITY / INSTITUTION NETWORK



### DATA FLOW

Afinion analyzers transfer data to RALS via a direct network connection, PC connection or a dataport, where it is evaluated and either sent directly to the EMR/EHR/LIS or held for review.

### DIRECT-TO-DEVICE INTERFACING

Afinion analyzers transfer a wide range of data to RALS, including:

#### Results

- Patient and Flagged

#### Sample Types

- Patient
- Level 1 & 2 QC Controls

#### Analytes

- HbA1c; ALB; Creat; ACR

### HOW THE DEVICE CONNECTS

*Afinion 2:*

- Direct network connection

*Afinion AS100:*

- Direct network connection
- RALS Dataport
- RALS Remote Connect (RRC)

### BI-DIRECTIONAL COMMUNICATION

*From device to RALS:*

- Patient and QC Results
- Patient ID and Operator ID
- Date/Time
- Location

*From RALS to device:*

- Operator ID
- Reagents/Tests Materials
- Device Configuration
- Date/Time

### AVAILABLE REPORTS FROM THE RALS SYSTEM

- Anomalous Results by Device
- Anomalous Results by Operator
- Device Statistics Summary
- Institution Statistics
- Levey-Jennings
- LIS Statistics by Location
- LIS Statistics by Operator
- Location Statistics
- Operator Pending Recertification Date
- Operator Pending Recertification Range
- Operator Statistics
- Operators with No Certifications
- Patient Results by Operator
- Patient Results by Patient
- QC Test Results
- RALS Device Log
- Test Materials

TO LEARN MORE, VISIT [RALS.COM](http://RALS.COM)

© 2019 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners. Any photos displayed are for illustrative purposes only. Any person depicted in such photos is a model. INFO247653 AI-COLL-0038 Rev. 5 06/19.

