



DEVICE CONNECTIVITY FOR
POINT-OF-CARE TESTING

RALS™ DEVICE INTERFACING

The NOVA Biomedical Stat Profile® Prime Plus® is a comprehensive, whole blood critical care analyzer that combines blood gases, electrolytes, metabolites, CO-Oximetry, and 32 calculated results in a simple, compact analyzer.

CONNECTIVITY TO RALS ENABLES USERS TO:


- Centralize all test results into the RALS system
- Transfer results to an EMR/EHR/LIS
- Email results to operators and managers
- Access and manage mis-identified patient and out of range QC results



VISIT [RALS.COM](https://www.rals.com)

RALS™ | Nova Stat Profile® Prime Plus® Interfacing

HOSPITAL / FACILITY / INSTITUTION NETWORK



DATA FLOW

The Nova Stat Profile® Prime Plus® connects to RALS via NovaNet middleware.

COMMUNICATION

Bi-directional Communication:

Data sent from the device to RALS

- Patient Results
- Control Results
- Proficiency Results

Data sent from RALS to device via NovaNet

- Operator List

AVAILABLE REPORTS FROM THE RALS SYSTEM

- RALS Devices Log
- Workload Summary
- Anomalous Results by Device
- Anomalous Results by Operator
- Operator Pending Recertification Range
- Operator Pending Recertification Date
- Operators with No Certifications
- Operator Statistics
- Patient Results by Operator
- Patient Results by Patient
- LIS Statistics by Location
- LIS Statistics by Operator
- Institution Statistics
- Location Statistics
- Levey-Jennings
- QC Test Results
- Device Statistics Summary
- Test Materials

TO LEARN MORE, VISIT RALS.COM

© 2020 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners. Any photos displayed are for illustrative purposes only. Any person depicted in such photos is a model. AI-COLL-0211 Rev. 1 04/20

