





## New UI FAQ #3

### How do I add a New Location?

#### To add a New Location:

Navigate to the Locations page by clicking the Locations icon  from the menu on the left and click the Add  button.

#### 1. Enter information into the available fields

The two required fields\* are Institution and Name:

1

General

Institution \*

Name \*

Address

Phone

CLIA License

CAP License

LIS Memo

- The Institution is chosen from the list provided in the drop-down
- The Name is a free text field and should indicate an actual location within the institution such as ED.

Review the available optional fields, filling out anything necessary for your new location



## New UI FAQ #3

### How do I add a New Location?

## 2. General section:

### a. Aliases

2

Aliases

Facility Code

LIS Name

LDAP Alias

BACKNEXT

### b. ADT

3

ADT Mapping

Q Parent Location

Q Child Locations

Q ADT Units

BACKNEXT



## New UI FAQ #3

### How do I add a New Location?

#### c. Personnel

4 Personnel

Q Supervisor

Medical Director

Clinician

BACK NEXT

#### d. Notes

5 Notes

Notes

BACK SAVE




## New UI FAQ #3

### How do I add a New Location?

#### 3. Click Save or

After creating the new Location, remember to assign the following items to the new Location, if applicable:


1. **Operators**- Operators will need the new location in their list of device access locations, **if** the device interfaced to RALS can receive an operator list
2. **Devices**-assign meters to the new location
3. **Materials** (if applicable) if materials are uploaded from RALS to the device, the location must be added to the appropriate materials
4. **Recertification Sets** (if applicable)-if recertification sets are in use, the new location must be assigned to the appropriate set



Click  for the following settings

5. **Device Configuration** (if applicable)-If the meter assigned to that location receives a configuration from RALS, assign that location to the appropriate configuration
6. **Evaluation Criteria** add the new location to applicable set
7. **Flag Configuration** add the new location to applicable set
8. **Patient Ranges** (if applicable) add the new location to applicable set
9. **Respiratory Configuration** (if applicable) add the new location to applicable set
10. **Comments** (if applicable) If comments are in use, add the new location to applicable set

For assistance with the following:

**Download Points** or **Wireless Configuration** please access the Help File  in the upper corner of RALS or call RALS TOLL-FREE SUPPORT LINE: 877.627.7257