



DEVICE CONNECTIVITY FOR
POINT-OF-CARE TESTING

RALS™ DEVICE INTERFACING

The Medonic™ M-series are fully automatic hematology analyzers intended for in vitro diagnostic testing of human blood samples under laboratory conditions. The Medonic M-series device module will support the standard features of other devices supported by the RALS System.


*Connectivity by the RALS System
enables users to:*

- Centralize all POC test results into the RALS system
- Transfer results to an EMR/EHR/LIS
- Email results to operators and managers
- Access and manage misidentified patient and out of range QC results
- Monitor and manage operator certifications



VISIT RALS.COM

RALS™ | Medonic™ M-series Connectivity



DATA FLOW

Medonic M-series analyzers transfer data to the RALS System using RALS Remote Connect (RRC) Software where it is evaluated and either sent directly to the EMR/EHR/LIS or held for review.

UNI-DIRECTIONAL COMMUNICATION

From analyzers to RALS System:

- Patient tests (CBC w/3-part diff)
- QC test results
- Calibration test results
- Operator ID
- Date/Time
- Comments
- Materials
- Specimen Types: Venous and Capillary Blood Samples

HOW THE DEVICE CONNECTS

- RALS Remote Connect

AVAILABLE REPORTS

- Anomalous Results by Device
- Anomalous Results by Operator
- Institution Statistics
- LIS Statistics by Location
- LIS Statistics by Operator
- Operator ID Compliance by Location
- Operator Pending Recertification Date
- Operator Statistics
- Operator Proficiency
- Patient ID Compliance by Operator
- Patient ID Compliance by Location
- Patient Results by Location
- Patient Results by Operator
- Patient Results by Patient
- QC Test Results
- Device Statistics Summary
- Test Materials
- Workload Summary

TO LEARN MORE, VISIT RALS.COM

© 2019 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners. Any photos displayed are for illustrative purposes only. Any person depicted in such photos is a model. INFO217701 AI-COLL-0066 Rev 2 06/19

